



# **Shaking Incubators**

Antech unique technology on the motor driving and tray balance controlled by the PID microprocessor enable the shaking more smooth, stable and safe.

Scope of application:Cell culture/Extraction experiment/Sample mixing

Dyeing and decoloring/Hybridization experiment

Solubility test/Plasmid purification/Bacterial suspension preparation

Elution process/Expressed protein



#### **Features**

- Vertical double-layer structure, small footprint and large bottle capacity.
- The original single-axis balanced drive technology is ultra-low and quiet.
- Intelligent acousto-optic alarm environment scanning micro-processing controller.
- The large LCD backlit LCD display shows the set parameters and measured parameters.
- The operating parameters are encrypted and locked to avoid human misoperation.
- Run parameter memory function to avoid tedious operation.



#### **General Control**

P.I.D microprocessor with LCD display shows icons and logical menus. Automatic run is available with non-volatile set point memory after power interruption.

### **Shaking Control**

Speed range: 30 to 300 rpm. Stroke: 25mm. Timer range: 0 to 999 hours 59 min.

Audible/Visible alarm for shaking speed deviation.

Automatic shaking speed adjustment in case of excessive workload such as unbalanced load placement, unusual vibrations caused from unstable floor or external shock.

## **Temperature Control**

Ambinet.  $+5^{\circ}$ C to 60 °C (for ATS-1102C/1102) Min. +4) to 60°C (for ATS-2102C/2102/1112B/2112B/2112D/2112F) Microprocessor PID control and auto-tuning. Temperature adjustment is 0.1°C, Audible and visible alarm for high and low temperature

#### **Accessories**

UV germicidal light keeps the chamber contamination free. LED light for users easy to observe shaking condition  $\,$ 



## **Flask Capacity**

Flask Capacity	50ml	100ml	250ml	500ml	750ml	1000ml	2000ml	3000ml	5000ml
ATS-1102C/2102C	48	48	30	24		12	_		_
ATS-1102/2102	104	104	56	44	_	30	_	_	_
ATS-1112B/2112B	164	164	90	70	_	36	30	_	_
/2112D/2112F									



## **Specification**

Product model	ATS-1102C	ATS-2102C	ATS-1102	ATS-2102	ATS-1112B	ATS-2112B	ATS-2112D	ATS-2112F	
Chamber volume	172L	172L	336L	336L	639L	639	639L	639L	
Shaking mode	Orbital	Orbital	Orbital	Orbital	Orbital	Orbital	Orbital	Reciprocal	
Shaking Speed	30~300RPM	30~300RPM	30~300RPN	30~300RPM	30~300RPM	30~300RPM	30~300RPM	30~300RPM	
Stoke	25mm	25mm	25mm	25mm	35mm	35mm	50mm	25mm	
Temperature Range	+5 °C ~60 °C	4~60℃	+5 °C ~60 °C	4~60 °C	4~60 °C				
Performance	Heated and Refrigerated Heated and Refrigerated Heated and Refrigerated								
Tray size (W*D)	500x350mm		734*	164mm	920*520mm				
Control	P. I. D Microprocessor								
Control mode	Fixed value or program (upto 8 segment)								
Display	LCD display								
Air Convection	Forced Air								
Ambient temperature	10-35 °C								
Temperature Accuracy	0.1 °C								
Temperature uniformity	≤± 1 ℃ (at 37 ℃)								
Stroke Accuracy	±1RPM								
Timer	0-500hours								
Tray Quantity	2								
Net weight	150KGS	165KGS	200	235KGS		45	6KGS		
Interior size(W*D*H)	600*400*650mm		850×520×760mm		1100*700*830mm				
Exterior size (W*D*H)	700*600*1250mm		950×740×1400mm		1200×850×1550mm				
Accessory	LED light/UV lamp								
Power	600W 1120								
Voltage	AC 200~240V 50H or AC 200~240V 60Hz								
Interior material	High-quality 304 stainless steel Mirror type								
Exterior material	High-quality cold rolled steel with electrostatic powder spraying								
Safety	"Audible and visible: High/low temperature, high/low shaking speed,								
	Electrical overload protection, Electrical leakge protection, Compressor overload protection"								

## **Order information**

FT001	T clamp,S/S,with spring retainer for 50mL flask	FT007	T clamp,S/S,with spring retainer for 750ml flask
FT002	T clamp,S/S,with spring retainer for 100mL flask	FT008	T clamp,S/S,with spring retainer for 1000ml flask
FT003	T clamp,S/S,with spring retainer for 150ml flask	FT009	T clamp,S/S,with spring retainer for2000ml flask
FT004	T clamp,S/S,with spring retainer for 250ml flask	FT010	T clamp,S/S,with spring retainer for 3000ml flask
FT005	T clamp,S/S,with spring retainer for 300ml flask	FT011	T clamp,S/S,with spring retainer for 5000ml flask
FT006	T clamp,S/S,with spring retainer for 500ml flask		