

Certificate of Compliance



No. 0P170328.QAS0017

Certificate's
Holder:

Qingdao Antech Scientific Co., Ltd.
No. 2 Yuantou Road, Qingdao, Shandong, China

Certification ECM
Mark:



Product:
Model(s):

Liquid Nitrogen Container
(see the following annex I)

Verification to:

Standard:
EN 13458-1:2002, EN 13458-2:2002+AC:2006,
EN 13458-3:2003+A1:2005

related to CE Directive(s):
2014/68/EU (Pressure Equipment)
ANNEX III MODULE A1

Remark: The product(s) has been verified on a voluntary basis. The product(s) satisfies the requirements of the Certification Mark of ECM, in reference to the above listed Standard(s). The above Compliance Mark can be affixed on the product(s) accordingly to the ECM regulation about its release and its use. The regulation can be found at www.entecerma.it. This Certificate of Compliance can be checked for validity at www.entecerma.it

This verification doesn't imply assessment of the production of the product(s).

Additional information, clarification about the **CE** Marking:



We attest that a TCF for the **CE** Marking process is in place. Whereas the Manufacturer is Responsible to start the **CE Marking Certification Procedure** through an appointed Notified Body and the perform all the necessary activities, as required by the Directive and accepted by the Notified Body, before placing the **CE** Mark on the product(s).

Date of issue 28 March 2017

Expiry date 27 March 2022

Chief Manager
Marta Morina

Deputy Manager
Amanda Payne

Annex I

No. 0P170328.QAS0017



Model(s):

CryoMajor 2, CryoMajor 2/30, CryoMajor 2/35, CryoMajor 3/50, CryoMajor 6/50, CryoMajor 10/50, CryoMajor 10/80, CryoMajor 10/125, CryoMajor 13/50, CryoMajor 15/50, CryoMajor 15/80, CryoMajor 15/125, CryoMajor 20/50, CryoMajor 20/80, CryoMajor 20/125, CryoMajor 20B, CryoMajor 25/50, CryoMajor 25/80, CryoMajor 30/50, CryoMajor 30B, CryoMajor 30/80, CryoMajor 30B/80, CryoMajor 30/125, CryoMajor 35/50, CryoMajor 35B, CryoMajor 35/80, CryoMajor 35B/80, CryoMajor 35/125, CryoMajor 47/127, CryoMajor 47/216, CryoMajor 50B/50, CryoMajor 50B/80, CryoMajor 50B/125, CryoMajor 80, CryoMajor 100, CryoMajor 130, CryoMaster 75, CryoMaster 100, CryoMaster 125, CryoMaster 600, CryoMaster 750, CryoMaster 900, CryoMaster 2400, CryoMaster 3000, CryoMaster 3600, CryoMaster 4800, CryoMaster 6000, CryoCarrier 3/50, CryoCarrier 6/80, CryoCarrier 8/80, CryoCarrier 10/125, CryoCarrier 25/216, CryoCenter 5, CryoCenter 15, CryoCenter 30, CryoCenter 50, CryoCenter 100, CryoCenter 150, CryoCenter 175, CryoCenter 200, CryoCenter 240, CryoCenter 300, CryoCenter 500, CryoSmart 600, CryoSmart 750, CryoSmart 900, CryoSmart 2400, CryoSmart 3000, CryoSmart 3600, CryoSmart 4800, CryoSmart 6000, CryoMonitor 1000, CryoMonitor 2000, CryoMonitor 3000, CryoMansion 10K, CryoMansion 24K, CryoMansion 38K, CryoWideneck 10K, CryoWideneck 24K, CryoWideneck 42K, CryoWideneck 45K, CryoPatronus 22K, CryoPatronus 52K, CryoPatronus 100K, CryoPatronus 126K, CryoMatrix 13K, CryoMatrix 19K, CryoMatrix 26K, CryoMatrix 36K, CryoMatrix 43K, CryoMatrix 50K, CryoMatrix 59K, CryoMatrix 76K, CryoMatrix 88K, CryoMatrix 95K, CryoMatrix 121K, CryoMatrix 128K, CryoBiobank 22K, CryoBiobank 24K, CryoBiobank 52K, CryoBiobank 126K ZYB-5, ZYB-8, LNM-1000, LNM-2000